

UNICO® Service Parts List for Model 1200 and 1201 Spectrophotometer

Parts #	Description
S-1200-401	UNICO Application Software for PC's, Window's 95 or 98, 2000, Windows ME, Windows XP required. Programs include Standard Curve, Abs./%T/Conc., and Abs.vs.Time. Includes S-2200-226 RS-232 serial cable.
S-1200-101	Test Tube V-type Cell Holder with base for 8~25 mm diameter test tubes
S-1200-102	Single Square Cuvette Holder with base for 10 mm Cuvettes
S-1200-103	4-position Cell Holder for 10 mm square cuvettes (Standard in model 1201)
S-1200-104	Longpath 4-Cell holder for up to 50 mm rectangular cuvettes
S-1200-105	Water-jacketed cell holder with base and panel for temperature control
S-1200-106	Universal Base for V-type holder, single cell holder and water-jacketed holder
S-1200-111	Black Block to set 0%T
S-2100-115	Holmium Oxide Filter, 10mm x 10mm x 45mm
S-2100-116	Didymium Filter, 10mm x 10mm x 45mm
S-1100-206	RS-232 Printer (Requires S-1100-207 Cable)
S-1100-207	Printer Cable for RS-232 printer (Male 25-pin to Female 9-pin, Null modem)
S-1200-220	Power cord, European plug
S-2100-226	RS-232C Cable for software (9-pin to 9-pin, Female/Female, Null Modem)
S-90-301	Test Tube Cuvettes, 10 mm diameter, Box of 12
S-90-302P-100	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Box of 100
S-90-302P-500	Disposable Cuvettes, Polystyrene, 10 mm pathlength, Case of 500
S-90-304G	Square Cuvettes, Optical Glass, Set of 2
S-90-305P	Disposable Semi-Micro Cuvettes, Polystyrene, 10 mm pathlength, 500pcs.
S-1100-101	Experiment Manual includes safety in the lab, experiments, & instructor's guide
S-1200-505	Tungsten-Halogen Lamps, Package of 1 (6V/10W, G4 pin)
S-1100-508	Printer Paper, Package of 3
S-1200-511	User Manual
S-1200-512	Dust Cover
S-1200-513	Fuse, 3 A, Package of 1
S-1200-520	Service Manual
S-1200-505	Halogen Bulb (6V 10 W G4 type) Package of 2
S-1200-535	Fuse, 3 A, Quantity 1 (size 5 x 20)
S-1200-801	Power Cord USA Plug
S-1200-801E	Power Cord European Plug
S-1200-802	On/Off Switch
S-1200-803	Lamp Socket (G4 type)
S-1200-804	110V/220V Selection Switch
S-1200-806	Transformer
S-1200-811	Detector/Amplifier Board
S-1200-812	CPU Board with Digital Display
S-1200-813	Power Supply Board
S-1200-891	Wavelength Knob
S-1200-892	Rubber Feet
S-1200-893	Sample Compartment Front Panel
S-1200-894	Operation Panel
S-1200-895	Wavelength Indicator
S-1200-896	4-position Cuvette Holder Rod
S-1200-897	Wavelength Disc

United Products & Instruments, Inc.
 182 Ridge Road, Suite E
 Dayton, NJ 08810, USA
 Tel: 732-274-1155 Fax: 732-274-1151
<http://www.unicol.com>
 E-mail: sales@unicol.com

United Products & Instruments, Inc.
182 Ridge Road, Suite E
Dayton, NJ 08810, USA
Tel: 732-274-1155 Fax: 732-274-1151
<http://www.unicol.com>
E-mail: sales@unicol.com