



UNITED PRODUCTS & INSTRUMENTS, INC.

SpectroQuest 2800 UV/Vis Spectrophotometers Parts

Part Number	Accessories and Components
SQ2800-101P	Test Tube Holder Kit for 8–20 mm diameter test tubes. Includes Universal base, one V-type round tube holder
SQ2800-102P	Long Path Cell Holder Kit for 4 cells up to 100 mm. pathlength. Includes Universal Base & one holder
SQ2800-102P-50	Long Path Cell Holder Kit for 4 cells up to 50 mm. pathlength. Includes Universal Base & one holder
SQ2800-104P	Cylindrical Cell Holder Kit for cells up to 100 mm pathlength. (22mm dia.). Includes Universal Base & one holder
SQ2800-105P	Water-Jacketed Cuvette Holder Kit for 10 mm single square cuvette Includes universal base, water-jacketed holder with flowthru panel (Water Bath Required)
SQ2800-106P	Micro Cell Holder Kit for use with 10 mm microcuvettes. Includes adjustable x-y base and single cuvette holder for 10mm square cuvette
SQ2800-107	Peltier controller unit for controlling temperature from 15°C to 40°C (requires S/N2800-107A peltier thermoelectric cell holder with panel)
SQ2800-107A	Peltier Thermoelectric Single Cell Holder with panel
SQ2800-107P	Peltier Package for continuous temperature control from 15°C to 40° C. Includes p/n SQ2800-107A peltier thermal electric 10mm cell holder with panel and p/n SQ2800-107 peltier controller unit
SQ2800-108A	Flow-thru Panel for use with ambient flowcell system (Requires flowcell, tubing and p/n SQ2800-109 flow-thru/peltier controller)
SQ2800-108P	Ambient Flow-thru Package. Includes p/n SQ2800-109 Flow-thru/Peltier Controller and p/n SQ2800-108A flow-thru panel (requires flowcell and tubing)
SQ2800-109	Peltier/Sipper Controller Only (Thermoelectric Temperature/Flow-thru Controller)
SQ2800-109A	Peltier/Sipper Panel and thermoelectric single cell holder for temperature-controlled flow-thru setup (Requires flowcell, tubing and p/n SQ2800-109 peltier/sipper controller)
SQ2800-109P	Peltier/Sipper Package for flow through and continuous temperature control from 15°C to 40°C Includes p/n SQ2800-109 peltier/flow-thru controller and p/n S-2800-109A peltier/sipper cell holder with panel (requires flowcell and tubing)
SQ2800-115	Holmium Oxide Filter
SQ2800-116	Didymium Filter
SQ2800-121	8-Position Auto Cell Changer for use with SQ2800/2802/3802 Series
SQ2800-131	Reflectance Measurement Attachment
SQ2800-132	Thin Film Measurement Attachment
SQ2800-226	Computer RS-232 Cable for use with application software (9-pin to 9-pin Female/Female, Null Modem)
SQ2800-401	Advanced PC Window software package featuring %T/Abs/Conc, Standard curve, Kinetics, DNA/Protein, Scanning, Multi-wavelength. RS232 serial cable and manual included. For SQ2800
SQ2800-405	Data Download-to-PC software for SpectroQuest Standalone models SQ2800/SQ2802 Kinetics, DNA/Protein, Scanning, Multi-wavelength. RS232 serial cable and manual included.
SQ2800-515	Tungsten-Halogen Lamp for SQ2800 Only. package of 1
SQ2800-525	Deuterium Lamp with bracket for SQ2800 series only, package of 1
SQ2100-535	Fuse, 3 A, Quantity 1