

Part List of PowerSpin M24

Parts No.	Specification	Qty Req.
C800-02	Regular tube shield	24
C800-03	Adapter for short tubes 4-5ml, pack of 10	
C800-04	Pediatric and microtube adapter, pack of 2	
C-M24-01	24 place rotor	
C800-10	Fuse, 250V 4A	
C-M24-101A	Lid assembly*	(1set)
C-M24-101-1	Lid, molded plastic	1
C-M24-101-2	Pins, mounting lid to case	2
C-M24-101-3	End cap for mounting pins	2
C-M24-101-4	Clear window on lid, plastic	1
C-M24-101-5	Lock Pin on Lid	1
C-M24-101-6	Nut to secure lock pin on lid	1
C-M24-102A	M24 Housing assembly* complete with panel, display, on/off switch and lock	(1set)
C-M24-102	Electric Fan	1
C-M24-103	Switch, On/Off, black	1
C-M24-104	Inlet power socket with fuse holder	1
C-M24-105	Control Panel	1
C-M24-106	Display Module	1
C-M24-107	Power Indicator (Red LED)	1
C800-108A	Lock assembly*	(1set)
C-M24-108-1	Motor, lock	1
C-M24-108-2	Micro switch, motor positioning(lock)	1
C-M24-108-3	Micro switch, safety (lock)	1
C-M24-109	Nut and washer to secure the rotor	1
C-M24-110	Guard Bowl, stainless steel	1
C-M24-111	Rubber gasket	1
C-M24-112	Main Control Board (CPU board)	1
C-M24-113	Capacitor	1
C-M24-114	Brake Resistor	1
C-M24-115	Speed Sensor and Board	1
C-M24-115-1	Speed Sensor Cable and connector	1
C-M24-115-2	Speed Disc with black and white line	1
C-M24-116	Motor, AC 110V, 2P, 100W, 1.2A	1
C-M24-117	Motor support, upper cushion	3
C-M24-118	Motor support, lower cushion	3
C-M24-119	Motor support, metal frame	1
C-M24-120-1	Screws, 50mm long (motor support)	8
C-M24-120-2	Spring Washer, heavy duty (motor support)	8
C-M24-120-3	Nut (motor support, rubber feet)	12
C-M24-121	Metal base	1
C-M24-122	Rubber feet	4
C-M24-123	Screw and washer, housing-base mounting	12
C-M24-124	Power cord, 115V, USA	1
C-M24-125	Power cord, 230V, European	1

* Assembly represents a group of componensts for service ordering convenience.
It is NOT part of material list

United Products & Instruments, Inc.
 182 Ridge Road, Suite E, Dayton, NJ 08810
 Tel: 732-274-1155 Fax: 732-274-1151
<http://www.unicol.com>
[E-mail: sales@unicol.com](mailto:sales@unicol.com)