

UNICO® PowerSpin™ HX Horizontal Centrifuge

UNICO® PowerSpin™ Series tabletop centrifuges are economical, reliable and super quiet.

The new horizontal fixed speed PowerSpin™ HX spins at 3400 RPM. The horizontal rotor insures a clean, complete separation of cellular components yielding a superior serum or platelet poor plasma (PPP). Using a horizontal rotor for spinning Serum Separator Tubes (SSTs) creates a tight uniform horizontal seal, reducing the risk of red cell leakage and possibility of specimen deterioration. This is especially important when the specimen is sent out for testing.

The PowerSpin HX model C822 features a Zero RPM safety lock, so the lid cannot be open while the rotor is in motion. Digital timer and LCD display provide user with the precise control of the spin cycle.

Powerful 100W brushless motor is maintenance free for years of dependable use. Stainless steel guard bowl adds even more protection, and also is easy to clean and maintain.

Please contact us for additional information



Features and Specifications:

- Horizontal Six place rotor for 2 to 10 ml tubes
- With adapters tubes accept 2 ml, or 4-5 ml tubes
- Quiet at full speed, only 58 decibels
- Stainless Steel rotor chamber is easily accessible for cleaning and disinfection
- HX model runs at fixed 3400 RPM speed
- 0-99 minute digital timer (Model C822)
- Zero RPM Safety Lock (Model C822)
- Suction cup rubber feet add stability while reducing vibration
- On/off switch located on the rear panel
- Detachable power cord with fuse holder
- Exterior is chemical resistant and easy to clean and disinfect
- Lid window for strobe tachometer RPM check
- Latching lid with safety switch triggered to shut-off power to the motor when the lid is raised
- Standard 115V/60Hz

Ordering Information:

UNICO Model	Item #	Speed (RPM)	Timer	Max. RCF	Capacity	Horizontal Rotor	Motor	Max. Tube Size
PowerSpin™ HX	C822	3400	Digital	1680 g	6 x 10ml	6 place	AC brushless	16 x 110mm

Accessories and parts:

C822-01	6 place swing out rotor	Centrifuge Size:	14.4 in x 11.6 in x 8.9 in (36cm x 29cm x 22 cm)
C822-02	Regular tube shield, pack of 2	Shipping Size:	16.5 in x 13 in x 10 in (41x 32x 25 cm)
C800-03	Adapters for 5 ml tubes, pack of 10	Net Weight:	22 lbs (10 Kgs.)
C800-04	Pediatric or 2 ml tube adapters, pack of 2	Shipping Weight:	25 lbs (11 Kgs.)
C800-10	Fuse 250A 4A, each	Warranty:	One year, parts and labor

United Products & Instruments, Inc.

182 Ridge Road, Suite E, Dayton, New Jersey 08810

© 2007 UNICO®

Tel: 732-274-1155 or 1-800-588-9776 Fax: 732-274-1151

E-mail: info@unicosci.com Web: www.unicosci.com