

T-Racks™ Tube Racks

The top shelf of each T-Rack has cross-referencing Alpha – Numeric designations so you can track the exact location of specific tubes in the rack.



The only rack you will ever need! Starting with a single 50 place T-Rack, you can attach as many racks together as you have room for! These patented racks quickly and easily attach together with their unique built-in T-Connector, which allows you the flexibility to attach one, two, three, or however many racks you decide together!

How do I attach T-Racks together?

It's so easy, even a hospital administrator could do it!

T-Rack's T-Connectors have a positive "click-stop", which locks the racks together, so they will not slide apart while being handled, or moved around the bench top, or refrigerator, etc. This click-stop holds them tightly together, yet is very easy to disassemble as well!

To disassemble, simply slide the racks apart and use as single rack, or join together with other T-Racks.

You can attach 10-13 mm racks to 13-16 mm racks! Attach any configuration to fit your needs! T-Racks allow you the flexibility to select how best to use them in your lab.

T-Racks can be used to help with Accessioning and are a great batching aide.

Item No.	Description	UOM
54500	T-Rack 50 Place, 13-16 mm, Expandable	Each
54500-2	T-Rack 50 Place, 13-16 mm, Expandable	Pk/2
54510	T-Rack 50 Place, 10-13 mm, Expandable	Each
54510-2	T-Rack 50 Place, 10-13 mm, Expandable	Pk/2
54520	T-Rack 50 Place, (25 holes of 13-16 mm & 25 of 10-13 mm), Expandable	Each
54520-2	T-Rack 50 Places, (25 holes of 13-16 mm & 25 of 10-13 mm), Expandable	Pk/2

Storage Tray



Organize and store supplies, tube racks (as shown) and more in this three-section, see-thru tray. Great for batching samples and organizing your work flow, or to store materials in the refrigerator, etc. Made of durable 3/16" acrylic.

Item No.	Description	UOM	Size
52030	Storage Tray, 3 Section	EA	12.75"W x 10"D x 3"H



Four T-Racks shown joined together, both 10-13mm and 13-16mm racks are attached. Mix and match to meet your needs, or to help with Accessioning.

Specifications:

- Each T-Rack holds up to 50 tubes of varying sizes in two 25 hole sections
- The top shelf of each T-Rack has cross-referencing Alpha – Numeric designations so you can track the exact location of specific tubes in the rack
- Attach additional T-Racks together, make a single large rack, as big as you need!
- Made of durable, autoclavable polypropylene
- T-Racks will not break, or become brittle in low temperature refrigerators
- Easy to clean and disinfect
- Dimensions: Each 50 Place T-Rack measures: 8.75" x 3.75" x 5.25"
- Your T-Racks will arrive disassembled and snap them together in seconds

Each T-Rack holds up to 50 tubes of varying sizes in two 25 hole sections.

