



Combination Glove Box and Paper Towel Dispenser

Two-in-One! This versatile dispenser holds gloves and paper towels in one convenient place. We recommend gloves on top and towels on the bottom (to prevent dripping on the gloves with wet hands), but it's your choice - either slot can be used for gloves or towels. Made of sturdy 3/16" acrylic.

Item No.	Description	Size
22650	Glove/Paper Towel Dispenser	11" L x 4" D x 16.75" H

Culture Plate Rack and Caddy



Organize, batch, transport and incubate, all with one rack! our Stainless Steel wire racks are stackable and autoclavable. Each rack holds up to 48 culture plates. Stainless steel rack caddy safely transports two Plate Racks side by side and is also stackable for more efficient use of incubator space.

Item No.	Description	Size
44510	Culture Plate Rack	12.5"W x 5"D x 9.5"H
44525	Culture Plate Rack Caddy	13"W x 9.5"D x 11.25"H

4-Drawer Organizers

48510



48512

48513

Organize your work with these handy desktop units. Select 8" or 12" width, each with 4 drawers so you can see-through. Each drawer features movable dividers to customize the space to fit your needs. Great for storing controls, calibrators or reagents in and out of the refrigerator.

- 4 drawers to store and organize your supply
- Movable divider in every drawer for custom fit
- Optional tray to fit 8" and 12" organizers together
- Constructed of durable 3/16" acrylic

Item No.	Description	Size
48510	Organizer, 4-drawer, 8" wide	8"W x 5.25"D x 8.5"H
48512	Organizer, 4-drawer, 12" wide	12"W x 5.25"D x 8.5"H
48513	Tray to hold 48510 and 48512	20"W x 5.25"D x 1.5"H
48514	Combo of 8", 12" 4-drawer and tray	20"W x 5.25"D x 8.5"H

Kimwipes® Dispensers



The versatile large dispenser holds large or extra large Kimwipes® boxes, or 3 small boxes lined up together.

- Free up countertop space by wall-mounting Kimwipes
- Small and Large box dispensers available
- Includes mounting hardware
- Made of sturdy 3/16" acrylic

Item No.	Description	Size
20650	Dispenser, Small Box	5.25"W x 3"D x 4.75"H
20655	Dispenser, Large Box	15"W x 5.5"D x 4"H