

Tube Holders are "Vertical at Rest" to allow use of Cover Slips for Fecal Protocols.



PowerSpin™ Centrifuge, MX Horizontal Spin

UNICO is proud to present our powerful 6-place Horizontal Spin centrifuge, the PowerSpin™ MX! The PowerSpin™ MX spins at speeds of 1000 - 3400 RPM, in increments of 100 RPM. Horizontal swing-out rotor provides high G-force and flat separation meniscus line.

The all digital push button front panel makes it easy to select or modify speed and run time even while wearing gloves. The PowerSpin™ MX also features a two line digital display that shows status, time, speed (RPM), and RCF (g).

Each centrifuge is extremely quiet at 50 dB. The rubber shock absorbing suction cup feet hold the unit in place resulting in a quieter run with less vibration. Integrated True Zero RPM Safety System prevents lid from being opened while the rotor is in motion.

Specifications:

- 6 place swing-out horizontal spin rotor for 2 to 15 ml tubes
- Push button panel for convenience when wearing gloves
- Quiet running, only 50 decibels at 3400 RPM (full speed)
- Stainless steel rotor chamber is accessible for easy cleaning
- Easy to read digital display shows real time RPM and RCF readings
- Remembers last setting for quick and easy repeat testing
- Variable Speed, 1000-3400 RPM
- Suction cup rubber feet eliminate vibration
- 60 minute digital timer with tone signal
- Lid window for speed verification with strobe tachometer
- Zero RPM Safety Latching System to prevent lid from being opened while rotor is spinning
- Standard 115V/60Hz AC
- One Year Limited Warranty

C8606 Model MX Horizontal Centrifuge

Actual Size: 13.75" x 15" x 11"
Shipping Dimension: 18" x 19" x 15.5"
Net Weight: 31 lbs
Shipping Weight: 36 lbs

Item No.	Description
C8606	MX, 6 Place Horizontal Spin Centrifuge
C8606-01	6 Place Swing Out Rotor
C800-02	Tube Shield, Pkg. of 2
C800-03	Short Tube Adapters, for 4-7 ml tubes, Pk/10
C800-04	Small Tube adapters, for 2 ml tubes, Pk/2
C800-10	Fuse for PowerSpin Centrifuges

PowerSpin™ Centrifuge, HX Horizontal Spin

The horizontal fixed speed PowerSpin™ HX centrifuge spins at 3400 RPM. The horizontal rotor insures a clean, complete separation of cellular components. Using a horizontal rotor for spinning Serum Separator Tubes (SSTs) creates a tight uniform horizontal seal, reducing the risk of red cell leakage and possibility of specimen deterioration.

The C822 features a Zero RPM safety lock, so the lid can not be opened while the rotor is in motion. A digital timer and LCD display provide user with the precise control of the spin cycle. Powerful brushless motor is maintenance free for years of dependable use. Stainless steel guard bowl adds even more protection, and also is easy to clean and maintain.

Specifications:

- Horizontal Six place rotor for 2 to 10 ml tubes
- Quiet at full speed, only 58 decibels
- Stainless Steel chamber for easy cleaning and disinfection
- HX models run at fixed speed 3400 RPM
- 0-99 minute digital timer
- Zero RPM Safety Lid Lock
- Suction cup rubber feet add stability while reducing vibration
- Lid window for strobe tachometer RPM check
- One Year Limited Warranty

C822 Model HX Horizontal Centrifuge



Actual Size: 14.4" x 11.6" x 8.9"
Shipping Dim.: 16.5" x 13" x 10"
Net Weight: 22 lbs
Shipping Weight: 25 lbs

Item No.	Description
C822	HX, 6 Place Horizontal Spin Centrifuge
C822-01	6 Place Swing Out Rotor
C822-02	Tube Shield, Pkg. of 2
C800-03	Short Tube Adapters, for 4-7 ml tubes, Pk/10
C800-04	Small Tube adapters, for 2 ml tubes, Pk/2
C800-10	Fuse for PowerSpin Centrifuges