



Rock-IT Tube Rocker

Actual Size: 13.7" x 4.7" x 5.0"
 Shipping Dimension: 15.7" x 6.7" x 6.9"
 Net Weight: 4.8 lbs
 Shipping Weight: 5.8 lbs

Rock-IT™ Tube Mixers

UNICO models LTR100 and LTR200 Rock-IT Tube Rockers are heavy-duty, yet nearly silent mixers designed for years of everyday use. Smooth and gentle movement at 25 RPM assures a uniform, thoroughly mixed specimen, while preserving the integrity of the fragile blood cells. Get your hematology controls ready to use in no time for accurate results. The reversible non-slip pad is autoclavable and will accommodate larger and hard to mix tubes. Rock-IT's are an economical choice for your lab; let the Rock IT work for you!

- Capacity of 8 or 15 tubes, up to 120 mm long
- 25 rocking cycles per minute, for gentle mixing
- Non-slip, reversible and autoclavable tube pad
- Compact size uses minimal benchtop space
- One year limited warranty

Item No.	Description
LTR100	Rock-IT Jr. Tube Rocker, 8 Place
LTR100E	Rock-IT Jr., same as above, 220V
LTR200	Rock-IT Tube Rocker, 15 Place
LTR200E	Rock-IT, same as above, 220V
LTR100-1	Replacement Tube Pad, 8 Place
LTR200-1	Replacement Tube Pad, 15 Place



Rock-IT Jr. Tube Rocker

Actual Size: 8.6" x 4.7" x 5.0"
 Shipping Dimension: 10.6" x 6.7" x 6.9"
 Net Weight: 4.0 lbs
 Shipping Weight: 5.0 lbs

Platform Lab Rotator



The LRT28 rotator is a compact, durable, variable speed orbital mixer. Accommodates a variety of printed or glass well slides, culture plates, flasks and small test tube racks, as well as RPR and VDRL slides or 96 well plates.


Specifications:

- 11" x 11" mixing plate with removable rubber pad
- Variable speed 25-250 rpm, orbital stroke of 20 mm
- 30 minute timer
- 3 position switch: Constant on/off/timed mixing
- Analog tachometer (adjustable)
- Small flask hold down spring included
- One Year Limited Warranty

Item No.	Description
LRT28	Variable Speed Laboratory Rotator

Shake-It! Lab Pipette Shaker

This type of shaker has been a reliable standard in the laboratory for generations. The 6-position Shake-It! is ideal for use in hematology, for special tests and manual techniques.



Specifications:

- Room for 6 pipettes
- 15 minute timer
- Triple position switch: On/Off/Continuous Use
- Acid/reagent resistant paint finish
- 110V AC
- One Year Limited Warranty

Item No.	Description
LP600T	6-Position Pipette Shaker