

Alcohol Dispenser

One-handed operation with our Alcohol Dispensing Pump Bottle; dispense onto pads, cotton balls, gauze, or brushes. Push down pump lid to dispense alcohol; the swing open lid reduces evaporation and odor. Holds 8 oz. alcohol (alcohol not included). Made of HDPE, sold as a case of 12 bottles.



Item No.	Description	UOM	Size
28050	Alcohol Dispenser	CS/12	2.5"W x 2.5"D x 5.75"H

Control Freek™ Control Racks



**1-Tier Rack
28071**



**2-Tier Rack
28072**



**Single Section
28070**

A great help for organizing your control vials, select a ONE, or TWO tier swiveling rack! Each rack comes with four or eight pizza slice shaped removable wedges. The Control Freek allows you to organize controls and calibrator chemicals in a manageable way – each rack swivels like a Lazy Susan for quick access. Each pizza slice is individually colored and can be removed to sit next to a particular analyzer, or can be independently put into a refrigerator, etc. Made of durable, easy to clean acrylic.

Item No.	Description	UOM
28071	Control Freek, 1 Tier Rack	EA
28072	Control Freek, 2 Tier Rack	EA
28070	Control Freek, Single Section	EA



See page 10 for our selection of Tube-CUBE's

Flexi-Phleb Phlebotomy Tray

The Flexi-Phleb has large carrying capacity without bulk. The tray provides a handy pull out drawer that is great for storing needles, bandages and other items. The top level accepts removable Tube-CUBE's (see page 10) that secure a variety of tubes, ice cups for temperature sensitive samples, or sharps containers. Back at the lab, merely exchange your collected samples in Tube-CUBE's at Specimen Processing with empty ones and you are ready to go again! Tray is constructed of tough, durable LDPE for long service and easy cleanup.

Flexi-Phleb includes two 9 place, one 4 place, and one 1 place Tube-CUBE's (sharps and alcohol dispenser not included).

Item No.	Description	UOM	Size
48700	Phlebotomy Tray	EA	9"W x 6"D x 14"H

Differential Counters

UNICO offers the 3 most popular counters: 2 key, 5 key with total window, and 8 key with total window. The 5 and 8 key models are sold with a WBC maturation series picto-strip just above the key windows.

Features:

- Three models: LBC3 - 2 key; LBC6 - 5 key; LBC9 - 8 key
- Heavy duty, impact resistant plastic housing
- Dual reset knobs
- Wide base with rubber feet for stability
- Signal bell when 100 cells have been counted



Item No.	Description
LBC3	2 Key Differential Counter
LBC6	5 Key, One Total Window Differential Counter
LBC9	8 Key, One Total Window Differential Counter