



Model C-MH30 Micro-hematocrit Centrifuge



Unico Model C-MH30 is made of formed metal and has a baked paint application that is acid and reagent resistant. A safety switch built into the handle assembly disconnects the power to the motor whenever the latch is lifted. The hematocrit tray is made of plastic composite material backed by an aluminum shell for durability and lighter weight. The brushed motor is mounted on rubber mounts to provide quieter running and less vibration.

Features:

- Twenty-four (24) place rotor with replaceable rubber ring and screw down metal lid.
- Nominal speed is 12,000 RPM
- Dependable mechanical 15-minute timer with bell signal at end of run
- Air-cooled, brush type AC motor
- All metal housing for years of dependable use
- Safety switch in latch shuts off motor when lid is lifted
- Electric motor braking
- Indicator light tells when unit is in use

Size	11.6 in x 10.3 in x 9.3 in (295 mm x 261 mm x 235 mm)
Shipping Size	12.7 in x 12.7 in x 12.2 in (323 mm x 323 mm x 311 mm)
Weight	17.4 lbs (7.9 kg)
Shipping Weight	19.4 lbs (8.8 kg)
Limited Warranty	One year

United Products & Instruments, Inc. 182 Ridge Road, Suite E, Dayton, New Jersey, 08810, USA
Tel: 1-800-588-9776 Fax: 732-274-1151 E-mail: sales@unicosci.com Website: <http://www.unicosci.com>