

UNICO® Colpo-Master™ III Colposcope

Triple Magnification, LED Illuminated Center Post Colposcopes

Introducing the NEW CS-305 Colpo-Master III Triple Mag Colposcope!

Simple and easy to use, featuring brilliant, sharp optics and three magnification settings: **3.75X, 7.5X** and **15X**. Changing magnification is a simple rotation of the knob, while staying in focus!

The state of the art LED illumination provides bright, white light that will not burn out, or burn too hot!

Features:

- Three Magnification Settings: **3.75x, 7.5x** and **15x**
- Rheostat intensity controlled, bright **LED** illumination
- Paired 12.5X Widefield Eyepieces with rubber eye-guards
- 3-D Stereoscopic Viewing
- Dual Diopter, independent eye focusing controls
- Flip-down green filter, built into the head
- 300 mm Working Distance
- Extra-Wide 115 mm Interpupillary Distance Adjustment
- Center Post stand with rolling, locking casters
- Light weight, easy to maneuver
- Broad Range Universal Power Supply 90-240V
- 7 Year Warranty, one year on electrical system
- Net Weight 23lbs
- Shipping Weight 32lbs



CS-305

UNICO ITEM #	DESCRIPTION COLPO-MASTER™ III	U/M
CS-305	Colpo-Master III Binocular LED Colposcope	1 ea

United Products & Instruments, Inc
182 Ridge Road, Suite E Dayton NJ 08810
800-588-9776 * info@unicosci.com * www.unicosci.com

Colpo-Master™ is a trademark of United Products & Instruments, Inc, dba UNICO

UNICO®