

UNICO[®] Colpo-Master[™] Video Packages

Analog and digital camera options

Introducing simple and effective options for physicians to capture images during exams or procedures for patient files, histories, or EMR programs. Provides real time video for exams, or as a teaching / training aide.

CS205T-LED
with CS-8185


Digital cameras offered by UNICO will capture crisp, high resolution color images utilizing custom software, developed with the physician office and hospital setting in mind. Captured images can be transferred to any EMR system for effective patient record management. Option of 3.0, or 5.6 Mega Pixel cameras. Connect using Windows XP or Windows 7 software included with the camera with USB cord and C-mount adapter.



Part	Description	Resolution	Sensor	Output
CS-8183	Digital Camera	3.0 MP	CMOS	USB
CS-8185	Digital Camera	5.6 MP	CCD	USB

A High-Resolution S-Video Color camera can provide a physician with the option of performing procedures while looking at a hi-res monitor instead of through the Colposcope for less eyestrain. The S-Video camera's Analog video signal can be converted to digital using the optional Frame Grabber (CS-8023).

Part	Description	Resolution	Sensor	Input/Output
CS-8021	S-Video Camera	480TVL	CCD	S-Video
CS-8022	Monitor	480TVL	NA	S-Video, BNC
CS-8023	Frame Grabber	NA	NA	S-Video/USB



CS103T-LED
with CS-8021

For your convenience we arranged several inclusive packages below. Please contact UNICO for additional information, or to schedule a rep. demo.

UNICO ITEM#	DESCRIPTION COLPO-MASTER VIDEO PACKAGES	U/M
CS103T-DV	Swing arm digital video package, includes CS103T-LED + CS-8185	1 ea
CS205T-DV	Center Post digital video package, includes CS205T-LED + CS-8185	1 ea
CS103T-SV	Swing Arm S-Video package, includes CS103T-LED + CS-8021	1 ea
CS205T-SV	Center Post S-Video package, includes CS205T-LED + CS-8021	1 ea
CS103T-SVM	Swing Arm S-video package with Monitor, includes CS103T-LED + CS-8021 + CS-8022	1 ea
CS205T-SVM	Center Post S-video package with Monitor, includes CS205T-LED + CS-8021 + CS-8022	1 ea
CS103T-SDV	Swing Arm S-video packaged with digital image grabber and monitor, includes CS103T-LED + CS-8021 + CS-8022 + CS-8023	1 ea
CS205T-SDV	Center Post S-video packaged with digital image grabber and monitor, includes CS205T-LED + CS-8021 + CS-8022 + CS-8023	1 ea

United Products & Instruments, Inc
182 Ridge Road, Suite E Dayton NJ 08810
800-588-9776 * info@unicosci.com * www.unicosci.com

ColpoMaster[™] is a trademark of United Products & Instruments, Inc, dba UNICO

UNICO[®]