

# G358 Dual Head Advanced Microscope

UNICO® is proud to announce the latest addition to its line: the newly designed G358 Dual Head Microscope. G358 Series has all the features and quality customers have come to expect from UNICO®. G358 comes equipped with a Dual Sliding Binocular Head. Both heads equipped with dual diopter adjustments; Standard High Eye Point 10X Wide Field Eyepieces allow user to work with the microscope without taking off glasses or protective masks.

- Sliding Binocular head with dual diopter adjustments. Heads rotate 360°
- WF10X High Eye-Point Eyepieces
- 4 advanced DIN standard achromatic objectives 4x, 10x, 40x, 100x (oil).
- Abbe condenser N.A 1.25 with iris diaphragm, filter holder and blue filter
- Full built-in mechanical stage, with hardened ball bearing mechanism, low coaxial controls and spring loaded slide holder. The stage is 142 mm x 133 mm in size and has an X-Y travel of 80 mm x 40 mm with graduation of 0.1 mm
- Ruggedly constructed metal stand and base features low mounted coarse and fine focus adjustment knobs for ease of use. The acid and alkali resistant paint is designed to provide years of durable use
- High intensity quartz halogen bulb (6V/20W) with rheostat to control intensity
- Koehler illumination field diaphragm for superior contrast control
- 5-year warranty excludes bulbs/fuses
- 40x Phase accessories available
- Semi-Plan and Plan Achromat objectives are optional
- Dust cover and extra fuses included



UNICO #	Head	Eyepieces	Objectives	Stage	Condenser	Illumination	Size	Weight
G358	Dual Binocular	10X WF	DIN 4X 10X 40X 100X Oil Achromatic	Mechanical	Abbe NA 1.25	Halogen 6V / 20W	18"x13"x10"	17 Lbs

**UNICO®**