

G380 Microscope Accessories

Here is a selection of some of the most commonly used accessories for the G380 Series. We inventory an extensive line of additional parts and accessories not shown here due to space limitations. If you do not see what you want or need, please contact us for more information (also see page 22 for more listings).



**G380-8907
Metal Carrying Case**

Our Metal Carrying Case is custom-designed to fit your G380 Series microscope. With foam-lined protection, it even has extra storage space for your optional LED illuminator and Semi-Plan Phase Kit. Perfect for a mobile practitioner. Complete with handle, lock and two keys.

Item No.	Description
G380-7001	G380 Series LED Battery Powered Illuminator (includes 3 AA batteries)
G380-7002	G380 Series Mirror Attachment
B6-2308	50X Oil Plan Objective
G380-8031	40X Semi-Plan Phase Objective and Slide-In Annulus
G380-8906	Darkfield Stop for G380 Series
B6-6010	Blue Filter, 31.5mm
B6-6011	Green Filter, 31.5mm
B6-6012	Yellow Filter, 31.5 mm
B6-8901	Eyepiece Reticule, 0.1mm divisions
B6-8903	Stage Micrometer Slide, 1mm in 0.01mm divisions
B6-8905	Immersion Oil, 7 ml
B6-7006	Halogen Bulb, 6V 20 Watt
B6-9005	Fuse, 1.5A, 250V
G380-8907	Metal Carrying Case for G380 Series



**G380-7001
LED Illuminator**

The LED Illuminator quickly turns your G380 into a cordless portable battery-operated field microscope. Includes 3 AA batteries. Installing the LED Illuminator is so easy, you can do it within 30 seconds. Works with any G380 Series microscope.



**G380-7002
Mirror Attachment**

Attaches to any G380 Series microscope. Quick and easy to install.



Optical Cleaning Kit

This kit is designed to protect and maintain one of the most valued assets in your laboratory. It has everything you need to take care of your microscope optics.

Includes:

- Blow Brush used to remove dust from lens surfaces
- Multi-coated formula Lens Cleaner
- Premium grade Lens Cleaning Tissue
- Micro-fiber Cloth
- Place Mats for eyepieces
- Carrying Case
- Detailed instructions

Item No.	Description
B6-8910	Optical Cleaning Kit



**G380-8031
SP Phase Kit**

Quickly turn your G380 Series microscope into a Phase Contrast microscope with our 40x Semi-Plan Phase Kit. Includes 40x Semi-Plan Phase Objective and Slide-in Annulus.

Did you know that UNICO services and maintains microscopes? Contact us for more details.