

IP758 Infinity Dual Binocular Microscope



IP758 Series Specifications:

- Infinity Corrected Optics
- Flat Field Achromat or Plan (PL) Achromat Objectives
- 4x, 10x, 40x and 100x (Oil) Infinity Objectives
- 40x and 100x Objectives with retractable tip
- Two pair (4), 10x Widefield, High Eyepoint eyepieces
- Dual Binocular Head, Infinity Corrected
- Dual Diopter Adjustments (4), one for each eyetube
- 55 – 75 mm Interpupillary Distance Adjustment
- Coaxial Coarse and Fine Focus control
- Low-Position Coaxial Stage Controls
- X-Y Axes Stage movement 80mm x 50mm
- Built-in Stage Tension adjustment
- Quartz Halogen 6 Volt / 20 Watt Illumination
- Intensity Rheostat for lighting control
- Transmitted Koehler Illumination with Graduated Iris
- 1.25 NA Abbe Condenser with Graduated Iris
- Built-in filter holder on condenser (blue, green, or yellow)
- Reverse-Position Revolving (4) place nosepiece
- Over-sized Mechanical Stage: 180mm x 150mm
- Dual Power Supply, 110V and 220V, built-in

Complete with extra Bulb, Dust Cover, Blue Filter, and Manual

IP758 Infinity Achromat, Dual Binocular Microscope

IP758PL Infinity Plan Dual Binocular Microscope