

L-TTR200D London Bus™, Double Decker Tube Rocker



The NEW UNICO® model **L-TTR200D** London Bus, Double Decker Tube Rocker has a capacity of 30 tubes on two, 15 tube shelves; and is an almost silent heavy-duty unit designed for years of everyday use. Available in both 110V and 220V versions. With this unit you can prep samples for diff counts, CBC's, mix injectables, or even prep analyzer controls.

Features:

- Capacity of 30 tubes or mix of tubes and vials (two 15 tube shelves)
- 25 gentle rocking cycles per minute
- Non-slip reversible and autoclavable tube pad(s)
- On/Off switch
- Will handle tubes up to 120 mm in length (up to 15 ml tubes)

Bench top Size	13.7 in x 4.7 in x 7.5 in
Shipping Size	16 in x 8 in x 7 in
Net Weight	5.2 lbs
Shipping Weight	7 lbs
Limited Warranty	One year