

UNICO Rock-IT™ Rockers Now with RED “PAUSE” button!



NEW! The RED "PAUSE" button available exclusively on all UNICO tube mixers. This red top-mounted button allows you to quickly pause mixing while you add or remove tubes. This handy feature extends the life of your rocker as you no longer grind gears while removing tubes. Newly redesign side panels add stability to the unit while retaining the small footprint on your benchtop.

Features:

- ◆ Option of 8, 15 or 30 tube platforms
- ◆ 25 gentle rocking cycles per minute
- ◆ Non-slip reversible and autoclavable tube pad
- ◆ Handles tubes up to 120mm in length

UNICO Model	Item Number	Capacity	Speed	Voltage	UOM
Rock-IT Rocker	L-TTR100	8 tubes	25cpm	110V	each
Rock-IT Rocker	L-TTR100-E	8 tubes	25cpm	220V	each
Rock-IT Rocker	L-TTR200	15 tubes	25cpm	110V	each
Rock-IT Rocker	L-TTR200-E	15 tubes	25cpm	220V	each
London Bus DD	L-TTR200D	30 tubes	25cpm	110V	each
London Bus DD	L-TTR200D-E	30 tubes	25cpm	220V	each



United Products & Instruments, Inc
www.unicosci.com * info@unicosci.com

182 Ridge Road, Suite E Dayton, NJ 08810
 Tel: 732-274-1155 * Fax: 732-274-1151

Rock-IT™ is a Registered Trademark of United Products & Instruments, Inc.