



# Programmable Multi Mix Rotator MTR 22



MultiMix MTR-22 provides several motion types in one instrument substantially increasing the options for mixing the examined materials. Our Multi-Shaker and Multi-Rotator range can be used for variety of applications, including the following: embedding blood cell suspensions, tissue-culture specimens, washing of precipitates, homogenizing blood solid with plasma, blot-hybridization.

- Multi Mix is an intelligent rotator for soft & hard mixing of micro-tubes and vacutainers.
- Several options in 1 instrument:
  1. Rotation
  2. Reciprocal rotation
  3. Vortexing
- Microprocessor controlled time, shaking speed, turning angle and pause
- General timer for control of program duration
- Standard platform capacity - 22 tubes up to 15 mm diameter (1.5 ml, 2 ml, 15 ml tubes)

<b>Technical Specifications</b>	<b>MTR 22</b>	Turning angle vortex mode	1-5
Capacity	22 tubes 1.8ml ~ 15ml	Vortexing times	1-5
RPM Range	5-30	Maximum timer	12 hours
Max Rotation Timer	999 sec	Power	110 / 220V
Max Reciprocation Timer	999 sec	Dimensions	342x177x192 mm
Turning angle reciprocation mode	15-90 ( 15 increment)	Weight	3kg

**UNITED PRODUCTS & INSTRUMENTS, INC.**

182 Ridge Road Suite E Dayton, NJ 08810 \* Tel (800) 588-9776 / (732) 274 1155 \* Fax (732) 274-1151

Web: [www.unicosci.com](http://www.unicosci.com) \* email: [info@unicosci.com](mailto:info@unicosci.com)