

FLEXI-PHLEB



Simple

- Compact Design
- Large Capacity

Flexible

- Phlebotomy
- Blood gas
- IV starting

Reliable

- Durable construction
- Light weight



A simple and compact design for the busy phlebotomist. The Flexi-Phleb has large carrying capacity without bulk. The tray provides a handy pull out drawer that is great for storing needles, bandages and other items. The top level accepts removable Tube-CUBE™ racks that secure a variety of tubes, ice cups for temperature sensitive samples, sharps containers, or whatever your specific activities demand. Back at the lab, merely exchange your collected samples in Tube-CUBE's at Specimen Processing with empty ones and you are ready to go again! Tray is constructed of tough, durable LDPE for long service and easy cleanup.

Flexi-Phleb Includes
4 Tube-CUBE's:
 2- 9 place
 1- 4 place
 1- 1 place
2 Squeeze Squeezer:
 Red Heart Phlebotomy Assistants (pn#28080-1)

Catalog #	Description	Dimensions
48700	Phlebotomy Tray	9" x 6" x 14"

The Tube-CUBE™ is an innovative and unique way to organize your tray. Available in 3 sizes to accommodate variety of tubes and accessories. No more ice spilled all over your basket or large tubes that have no place on a old fashion tube rack. Light weight durable plastic you can load up and easily carry a full tray.

Additional Tube-CUBE™ racks sold separately		
Catalog#	Tube-CUBE Description	Layout
54901	9 place for 8-16 mm tubes	● ● ● ● ● ● ● ● ●
54401	4 place for 26-30 mm tubes	● ● ● ●
54101	Single place for up to 66mm tubes	●

UNITED PRODUCTS & INSTRUMENTS, Inc
 182 Ridge Road, Suite E Dayton, New Jersey 08810
 Tel: 800-588-9776 * Fax: 732-274-1151
 E-mail: info@unicosci.com * Web: www.unicosci.com