

**New From UNICO!**

This beginner Stereo Microscope is designed for curious minds. Developed to be a “First Stereo” for a hobbyist or for K-12 education; the ST11 incorporates features and optics often found in more expensive microscopes. Slip-clutch focusing mechanism, locked eyepieces - all part of “student-proof” design. Adjustable diopter and interpupillary distance allow the ST11 to accommodate every viewer.

Rugged and dependable rack and pinion mechanism provides smooth and dependable precision focusing. Reversible Black & White stage plate, dust cover and user manual are included with every microscope.

The microscope works with ambient light. This makes it an ideal choice not only for the classroom (no cords to mess with and safe for young learners); it is also a perfect instrument for field work. Light weight and compact design will help you set up a workshop or after- school program in minutes.

UNICO microscopes are built to last and come standard with a 5 year warranty.

**Specifications**

- ❑ Binocular head with vertical eye tubes
- ❑ Coated glass optics
- ❑ 3-Dimensional viewing
- ❑ Locked-on, 10X widefield eyepieces
- ❑ Diopter and interpupillary distance adjustment
- ❑ Rubber eyecups
- ❑ Objectives: Paired 2X (20X Power)
- ❑ Rack and Pinion focus adjustment with slip clutch
- ❑ Working Distance: 71mm (2.75")
- ❑ Large base, acid and reagent resistant with locked-on stage clips
- ❑ Black/White reversible stage plate
- ❑ Non-Skid Rubber Feet, metal frame
- ❑ Dust cover and User's manual
- ❑ 5 year limited warranty

Actual Size: 6½" x 4 1/2" x 10 ¼"  
Shipping Dimension: 8" x 5 ¾" x 13"  
Shipping Weight: 2.5 lbs