

Tube Rocker

Product Images



Short Description

UNICO TEST TUBE ROCKER
UP TO 15 TUBES/VIALS CAPACITY
EASY OPERATION TOP PAUSE BUTTON
FIXED SPEED OF 25 ROCKING CYCLES PER MINUTE
POWER INPUT 220 V

Description

UNICO Rock-IT's now feature a PAUSE button, a useful feature designed to allow the user the opportunity to pause the rocking motion while adding or removing tubes. The feature also extends the life of the rocker motor by eliminating the need to hold the platform to add/remove tubes.

UNICO Rock-IT Tube Rockers are heavy-duty, yet nearly silent mixers designed for years of everyday use. Smooth and gentle movement at 25 RPM assures a uniform, thoroughly mixed specimen, while preserving the integrity of the fragile blood cells.

Get your hematology controls ready to use in no time for accurate results. The reversible non-slip pad is autoclavable and will accommodate larger and hard to mix tubes. Rock-IT's are an economical choice for your lab; let the Rock IT work for you!

Additional Information

Specification	Capacity	15 tubes
	Speed	25cpm
	Voltage	220V
Size and Weight	Size	8.6 in x 4.7 in x 5.0 in
	Shipping Size	10.6 in x 6.7 in x 6.9 in
	Weight	4.0 lbs (1.8 kg)
	Shipping Weight	5.0 lbs (2.3 kg)