

## Compact Incubator

### Product Images



## Short Description

---

UNICO COMPACT INCUBATOR  
6L CHAMBER CAPACITY  
REMOVABLE SHELF  
TEMPERATURE AMBIENT TO 60°C (140°F)  
POWER INPUT 110 VOLT

## Description

---

The L-CU60 is a highly accurate incubator that is small in size but big in performance. Featuring a see-through door that makes a visual check on your cultures quick and easy. With a capacity of 10 standard size culture plates per shelf, the nominal capacity is 20 plates. The acrylic body is a great insulator and is easy to clean and maintain. With a small footprint, this unit is perfect for labs with limited benchtop space

## Additional Information

---

Specification	Chamber Capacity	6 Liters (0.21 cubic feet)
	Chamber Size:Internal Dimensions	6.8