

Compact Incubator

Product Images



Short Description

UNICO COMPACT INCUBATOR
6L CHAMBER CAPACITY
REMOVABLE SHELF
TEMPERATURE AMBIENT TO 60°C (140°F)
POWER INPUT 220 VOLT

Description

The L-CU60 is a highly accurate incubator that is small in size but big in performance. Featuring a see-through door that makes a visual check on your cultures quick and easy. With a capacity of 10 standard size culture plates per shelf, the nominal capacity is 20 plates. The acrylic body is a great insulator and is easy to clean and maintain. With a small footprint, this unit is perfect for labs with limited benchtop space

Additional Information

Specification	Chamber Capacity	6 Liters (0.21 cubic feet)
	Chamber Size:Internal Dimensions	6.8