

## Advanced Laboratory Microscope

### Product Images



## Short Description

---

SEIDENTOPF TRINOCULAR HEAD  
WF10X PAIRED EYEPIECES  
4X,10X,40XR,100XR DIN PLAN ACHROMAT OBJECTIVES  
BUILT-IN MECHANICAL STAGE WITH ABBE NA1.25 CONDENSER & IRIS  
HALOGEN ILLUMINATION

## Description

---

The IP730/750 UNICO is designed to provide crisp, clear imaging on a versatile and ergonomic platform. Available with a variety of optics and configuration choices. The IP750 series is designed with high-end Infinity optics.

## Additional Information

Header	Bodies
IP730-3004	SEIDENTOPF BINOCULAR BODY WITHOUT EYEPIECES FOR IP730 SERIES MICROSCOPE
IP730-3005	SEIDENTOPF TRINOCULAR BODY, 30 DEGREE INCLINED WITHOUT EYEPIECES FOR IP730 SERIES MICROSCOPE
IP750-3004	SEIDENTOPF BINOCULAR BODY WITHOUT EYEPIECES FOR IP50 SERIES MICROSCOPE
IP750-3005	SEIDENTOPF TRINOCULAR BODY, 30 DEGREE INCLINED WITHOUT EYEPIECES FOR IP750 SERIES MICROSCOPE

Header	Optics
IP750-1502	10X EXTRA WIDEFIELD EYEPIECE FN 22MM FOR IP730 AND IP750 MICROSCOPE
IP730-2101	4X DIN ACHROMAT OBJECTIVE, N.A. 0.10
IP730-2102	10X DIN ACHROMAT OBJECTIVE, N.A. 0.25
IP730-2103	40X DIN ACHROMAT OBJECTIVE, N.A. 0.65, RETRACTABLE FRONT LENS
IP730-2104	100X DIN ACHROMAT OBJECTIVE, N.A. 1.25, RETRACTABLE FRONT LENS
IP730-8030	PHASE SET W/ PLAN PHASE 10X-40X-100X OBJECTIVES, TURRET CONDENSER WITH "0" POSITION
IP750-2101	4X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2102	10X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2103	40X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2104	100X INFINITY ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2301	4X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2302	10X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2303	40X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2304	100X INFINITY PLAN ACHROMATIC OBJECTIVE FOR IP750 INFINITY MICROSCOPE
IP750-2422	INFINITY PLAN PHASE OBJECTIVE 10X FOR IP750 SERIES
IP750-2423	INFINITY PLAN PHASE OBJECTIVE 40X FOR IP750 SERIES
IP750-2424	INFINITY PLAN PHASE OBJECTIVE 100X FOR IP750 SERIES
IP750-8030	PHASE SET W/ INFINITY PLAN PHASE 10X-40X-100X OBJECTIVES, TURRET CONDENSER WITH "0" POSITION

Header	Miscellaneous Accessories
IP750-5001	MECHANICAL STAGE WITHOUT ABBE CONDENSER FOR IP730 AND IP750 MICROSCOPE
IP750-6002	PHASE CONTRAST CONDENSER, TURRET, NA1.25 WITH 10X-40X-100X PHASE RING AND "0" BRIGHT FIELD POSITION
IP750-6003	DARKFIELD CONDENSER FOR IP730 AND IP750 SERIES MICROSCOPE
IP730-8002	DUST COVER FOR IP730 AND 750 SERIES MICROSCOPE
IP750-7006	6V 20W HALOGEN BULB, 4MM PIN FOR IP730 AND IP750 MICROSCOPE
IP750-8901	HARD CARRYING CASE FOR IP730 AND IP750 SERIES MICROSCOPE

Accessories

Specification	Head	Trinocular
	Eyepiece	WF10X
	Objectives	4x 10x 40XR 100XR Plan
	Stage	MECH.
	Condenser	Abbe N.A.1.25
	Illuminator	6V/20W /Halogen
Size and Weight	Size	10 in x 8.5 in x 15.4 in
	Shipping Size	11 in x 14 in x 18 in
	Weight	18 lbs
	Shipping Weight	20 lbs